

GE0-807^{Q&As}

Genesys Certified Professional 8 - System Consultant, Outbound Voice

Pass Genesys GE0-807 Exam with 100% Guarantee

Free Download Real Questions & Answers **PDF** and **VCE** file from:

https://www.geekcert.com/ge0-807.html

100% Passing Guarantee 100% Money Back Assurance

Following Questions and Answers are all new published by Genesys Official Exam Center

Instant Download After Purchase

100% Money Back Guarantee

- 🔅 365 Days Free Update
- 800,000+ Satisfied Customers





QUESTION 1

When the internal registrar is enabled and SIP Server is also configured to store registration information in the configuration database, where, in CME, can you find the contact option specifying current location?

- A. corresponding DN
- B. corresponding Agent Login
- C. corresponding Voice over IP Service
- D. corresponding SIP Server
- E. It is not visible in CME

Correct Answer: A

QUESTION 2

SIP Server starts successfully. However, its status on SCS is "Service Unavailable" and SIP Server distributes EventLinkDisconnected to clients. Which of the following circumstances would be responsible for the status and EventLinkDisconnected message?

- A. Stream Manager is not registered with SIP Server
- B. Media Gateway is not registered with SIP Server
- C. SIP Server fails to open the SIP port
- D. The network is down

Correct Answer: B

QUESTION 3

SIP Server supports active out-of-service detection for which of the following types of ONs? (Choose 2 answers)

- A. VoIP Service (MCU, treatment, etc.)
- B. Extension
- C. Trunk
- D. Route Point

Correct Answer: AB

QUESTION 4



Which of the following is SIP Server deployment modes used with a Framework 8 solution? Choose 3 answers

- A. Stand Alone Mode
- B. Application Server Mode
- C. Multimedia Mode
- D. Customer Side Proxy Mode
- E. CTI Mode

Correct Answer: ABD

QUESTION 5

Which statement is true?

A. Recording is not possible from a routing strategy.

B. If busy treatments are played from a routing strategy while Customer is waiting for an available Agent, as soon as a target is available, the busy treatment is stopped.

Correct Answer: B

GE0-807 Practice Test GE0-807 Exam Questions

GE0-807 Braindumps